



SIMPLEE®

All you need. No more. But no less.



ALL YOU NEED.

NO MORE.

BUT NO LESS.



ONE ADHESIVE FOR ALL OCCASIONS

SIMPLEE® UNIVERSAL BOND

Self-etching universal bonding, light-curing. Suitable for all etching techniques. Very good adhesion even on moist dentine. Can be used with all dual, self or light curing materials. Can also be used as a primer for metal and ceramics. Repairs to composite and ceramics.

Indications: All classes of direct restorations with composite (classes I to V). Fixation of root posts, crowns, inlays / onlays, veneers.

5 ml Bottle
156 110

€ X



SIMPLEE® PROTHESENFINISH

Gloss spray for prostheses and plastics. Take care of one uniform satin finish with mint scent.

Spray on from about 4 - 5 cm distance and polish lightly with a soft cloth.

75 ml

093 996 1 Can € X



SIMPLEE® SCAN SPRAY

Classic matte scan spray with fine but dense powder film of 5 µm for daily laboratory use.

Dries quickly after application. by increased density also good for masking metal suitable for restorations. TiO2 free. Color: White.

200 ml

093 994 1 Can € X



SIMPLEE® MIXING PAD

For impression materials, cements and composites.

Made of abrasion-resistant paper specially coated on one side, blocked on both sides.

7 x 9,5 cm, 60 Sheets

300 105 1 Pckg. € X

12,5 x 18 cm, 25 Sheets

300 106 1 Pckg. € X



 **SIMPLEE®**
ETCHING GEL

Etch gel with 36 % phosphoric acid.
For the conditioning of dental enamel and dentine with optimal thixotropic flow characteristics. It is used for etching tooth substance within the framework of the selective enamel etching technique and total-etch technique
Colour: blue-green.

20 x 1,4 ml Syringe,
25 Application cannulas
156 115 1 Pckg. €X

Application cannula 22G,
angle 45°, blue,
10 mm long, 100 Pcs.
156 104 1 Pckg. €X

Application cannula 20G,
angled, black,
11 mm long, 100 Pcs.
100 112 1 Pckg. €X



NEW



NEW



NEW

SIMPLEE®
SIMPLEEJECT

Hypodermic 3 ml syringe irrigation syringe combined with blunt (notched) cannula. Sterile packaged, assembled ready-to-use.

Size	100 Pcs.
23G, 0,6 mm	500 727
27G, 0,4 mm	500 728
1 Pckg.	€ X

SIMPLEE®
ALGINATE

Dust-free alginate for a detailed impression. With long dimensional stability (100 hours). Peppermint flavour. Shelf life 5 years from the date of manufacture.

Available in three versions:

Simplee Alginate:
normal setting, mixing time: approx. 45 seconds
preparation until impression taking: 1 minute 45 seconds
setting time in the mouth: 30 seconds
total setting time: 1 minute 55 seconds
colour: blue.

Simplee Alginate quick:
fast setting, mixing time: approx. 45 seconds,
preparation until impression taking: 1 minute 30 seconds
setting time in the mouth: 30 seconds
total setting time: 1 minute 40 seconds
color: blue

Simplee Alginate color:
normal setting, mixing time: approx. 45 seconds
preparation until impression taking: 1 minute 45 seconds
setting time in the mouth: 30 seconds
total setting time: 1 minute 55 seconds
colour: purple to pink to white..

	Bag: 500 g	Bulkpack: 20 bags à 500 g
Alginate	156 201	156 200
Alginate quick	156 203	156 202
Alginate color	156 205	156 204
1 Pckg.	€ X	€ X

SIMPLEE®
SIMPLEEDRY

SimpleeDry absorbs 30 times its own weight. Complete coverage of parotid orifice and entire buccal mucosa. No loose fibers or other residues left in mouth. Flexible, accommodates to cheek movements. Comfortably stays in place, easily removed with water spray

Size	50 Pcs.
small	156 552
large	156 553
1 Pckg.	€ X

Note: Simplee product line prices are non-discountable and exclude our regular orders and special discounts.




 **SIMPLEE®**
FACE MASK IIR

3-ply face mask with elastic loops and nose clip

Filter performance >98% according to EN 14683:2019 Type IIR

Very comfortable to wear, breathable and liquid-resistant.

Colour	50 Pcs
black	700 131
magenta	700 132
indigo	700 133
fossil	700 134
1 Pckg.	€ X





200 PCS



SIMPLEE[®]
NITRILE GLOVE
200 PCS

Powder-free nitrile glove with good grip thanks to special finger-tip texturing. AQL 1,5, manufactured according to EN455, also approved as a barrier glove according to EN 374. Ambidextrous fit, color: blue and white.

Size	200 Pcs. blue	200 Stck. white
XS	156 570	156 574
S	156 571	156 575
M	156 572	156 576
L	156 573	156 577
1 Pckg.	€ X	€ X




 **SIMPLEE®**
BIB HOLDER

Disposable self adhesive bib holder with logo.
250 pcs/pack.

Comply with hygiene recommendations, minimizes
requirements for reprocessing.

250 Pcs.

156 528 1 Pckg. € X





 **SIMPLEE®**
PATIENT BIB

Moist barrier protection during treatment
Manufactured with 2-ply paper and PE-coating


33 x 45 cm, 500 Pcs:

white		
156 525	1 Carton	€ X
dark blue		
156 525	1 Carton	€ X
green		
156 527	1 Carton	€ X



COMING SOON




 **PAPER TOWEL**
1-LAYER

natural:
Zickzack-fold, 24 x 21 cm,
5 000 Pcs. (250 pcs bundles)
156 535 1 carton € X

green (recycled material):
Zickzack-fold, 23 x 24,8 cm,
4000 Pcs. (200 bundles)
156 538 1 Carton € X



 **PAPER TOWEL**
2-LAYER

natural (cell stoff):
Layer-fold, 23 x 33 cm,
2,560 Pcs.
156 537 1 Carton € X


natural (cell stoff):
Zickzack-fold, 23 x 24,8 cm,
3,750 Pcs. (250 Bundles)
159 553 1 Carton € X



 **SIMPLEE®**
PROBES

Hollow cylindrical handle.

- Fig. 9
156 610 1 pc. € X
- Fig. 23
156 611 1 pc. € X
- Fig. 8
156 612 1 pc. € X
- Fig. 17
156 613 1 pc. € X



Nr. 19



Nr. 20



Nr. 21



SIMPLEE® HEIDEMANN SPATULA

Hollow cylindrical handle.

Nr. 19, 2,0 mm
156 605 1 pc. € X

Nr. 20, 2,5 mm
156 606 1 pc. € X

Nr. 21, 3,0 mm
156 607 1 pc. € X

Fig. 1/2



Fig. 3/4



Fig. 5/6



Fig. 7/8



Fig. 9/10



Fig. 11/12



Fig. 13/14



Fig. 15/16



SIMPLEE® GRACEY CURETTES

Gracey curettes for sub-gingival scaling.

Fig. 1/2
156 615 1 pc. € X

Fig. 3/4
156 616 1 pc. € X

Fig. 5/6
156 617 1 pc. € X

Fig. 7/8
156 618 1 pc. € X

Fig. 9/10
156 619 1 pc. € X

Fig. 11/12
156 620 1 pc. € X

Fig. 13/14
156 621 1 pc. € X

Fig. 15/16
156 622 1 pc. € X

Fig. 204S



Fig. 204SD



Fig. H6/H7



SIMPLEE® SCALER

Scaler for removal of subgingival tartar and plaque.

204S
156 627 1 pc. € X

204SD
156 628 1 pc. € X

H6/H7
156 629 1 pc. € X

NEW

 **SIMPLEE®**
TWEEZERS

Tweezers with serrated tip.

Meriam 16 cm

156 600 1 pc. € X

College 15 cm

156 601 1 pc. € X



 **SIMPLEE®**
MOUTH MIRROR
HANDLES

Hexagonal, hollow


156 630 1 pc. € X

Octagonal, massive

156 631 1 pc. € X

Finger shaped, hollow

156 632 1 pc. € X



 **SIMPLEE®**
ASPIRATION-
SYRINGE

Cylinder cartridge syringe (1,8 ml) with
tilt function, 3 aspirations hooks

156 625 1 pc. € X






SIMPLEE®
SCALER ANSÄTZE
 ZEG connection made of corrosion-resistant, flexible steel.

for Dentsply Sirona* :

Tartar removal devices:



sr3l Universal
156 013 1 pc. € X



sr4l Universal fine
156 014 1 pc. € X



srPE1 Perio
156 037 1 pc. € X



sr2l Versatile
156 036 1 pc. € X



sr1l General
156 035 1 pc. € X



srPE2 Perio right-angled
156 038 1 pc. € X



srPE3 Perio left-angled
156 039 1 pc. € X

Torque wrench:



for DS*
156 025 1 pc. € X

for KaVo®:

**KaVo® Sonicflex 2000 N/
Lux 2003 LX:**



k5 Universal
156 007 1 Stck. € X



k6 Sickle
156 008 1 Stck. € X



k7 Perio
156 009 1 Stck. € X



k60 long
156 031 1 Stck. € X



k61 paro
156 033 1 Stck. € X



k62 paro
156 034 1 Stck. € X

KaVo® Piezolux/ Sonosoft:



k11 Universal
156 010 1 pc. € X



k9 Sickle
156 011 1 pc. € X



k10 Perio
156 012 1 pc. € X

KaVo® Sonicflex quick 2008:



k5A Universal
156 003 1 pc. € X



k6A Sickle
156 004 1 pc. € X



k7A Perio
156 005 1 pc. € X



k8A Perio extra long
156 006 1 pc. € X



k60A long
156 032 1 pc. € X

KaVo® Piezo LED:



k201
156 000 1 pc. € X



k202
156 001 1 pc. € X



k203
156 002 1 pc. € X

Torque wrench:



for KaVo Sonicflex*
156 024 1 pc. € X



for KaVo Sonicflex*
156 029 1 pc. € X



for KaVo PiezoLED/Piezo soft*
156 030 1 pc. € X

for NSK*:

Tartar removal devices:



nS1 Universal
156 015 1 pc. € X



nS2 Sichel
156 016 1 pc. € X



nS3 Perio
156 017 1 pc. € X



nG1
156 040 1 pc. € X



nG16
156 041 1 pc. € X



nG4
156 042 1 pc. € X



nG6
156 043 1 pc. € X



nG8
156 044 1 pc. € X



nS20
156 045 1 pc. € X

Torque wrench:



for NSK*
156 026 | 1 pc. | € X

for ACTEON*:

Tartar removal devices:



s1 Universal
156 021 1 pc. € X



s10X interdental/ approximal
156 022 1 pc. € X



s10P sub-/ supragingival
156 023 1 pc. € X



S10Z subgingival
156 046 1 pc. € X



S1S sub-/supra slim
156 047 1 pc. € X



sH3 anterior Perio
156 048 1 pc. € X

Torque wrench:



for ACTEON*
156 028 1 pc. € X

for EMS*:

Tartar removal devices:



eA Universal
156 018 1 pc. € X



eP Perio
156 019 1 pc. € X



ePS Perio fine
156 020 1 pc. € X

Torque wrench:



for EMS*
156 027 1 pc. € X

SIMPLEE® WATTEROLLEN

Cotton rolls made of 100% cotton.
Skin friendly, does not stick to oral mucosa.

Size	300 g
1	156 565
2	156 566
3	156 567
1 Pkg.	€ X



NEW



SIMPLEE® BIO DISPOSABLE TOOTHBRUSH

Disposable toothbrush with wheat based handle.
4 filament rows impregnated with toothpaste.
Individually hygienically packaged.

blue, 100 Pcs.
728 817 1 Pkg. **€ X**




SIMPLEE® DAPPEN GLASS

Smooth and jagged plastic dappen glasses.
The Bio dappen glasses are made from fossil-
free mtrl (tree oil).

Cannot be sterilized. Not a medical device..

Colour	smooth, 50 Pcs.	jagged, 50 Pcs.
white	411 501	411 506
yellow	411 502	–
blue	411 503	411 507
red	411 504	–
black	411 505	–
1 Pkg.	€ X	€ X
Farbe	BIO smooth, 50 Pcs.	BIO jagged, 50 Pcs.
white	411 508	411 509
1 Pkg.	€ X	€ X





SIMPLEE® PROPHY BRUSH

With nylon filament. Not for sterilisation.
RA shank.

Hardness	100 Pcs.
pink, soft	156 805
white hard	156 806
1 Pckg.	€ X

SIMPLEE® PROPHY CUPS

Latex-free disposable prophy cups in 2
degrees of hardness. RA shank.

Hardness	24 Pcs.
blue, soft	156 810
yellow hard	156 811
1 Pckg.	€ X

SIMPLEE® WIPES RSD 400

Ready-to-use, alcohol-reduced disinfectant wipes
(18x20 cm) for the rapid disinfection of non-invasive
medical devices and treatment areas.
Very good material compatibility.
In a practical flow pack that can be resealed.
Broad spectrum of activity and rapid effectiveness.
Cloth material made from renewable raw materials.
Neutral scent

6 Flow-packs 100 wipes/pack
156 706 1 Pckg. **€ X**

SIMPLEE®
BRUSH INSERTS

100 Pcs.:

fine, transparent shaft

300 100 1 Pckg. € X

medium, red shaft

300 101 1 Pckg. € X

large, transparent shaft

156 103 1 Pckg. € X

Holder for Simplee Micro Tips/

Brush inserts

156 101 1 Pckg. € X



SIMPLEE®
MICRO TIPS

Miniature application tips, flocked.
For the application of liquids and gels.

transparent, 100 Pcs.

156 100 1 Pckg. € X

Holder for Simplee Micro Tips/

Brush inserts

156 101 1 pc. € X



SIMPLEE®
FLEXIBRUSH

Bendable brush for application
of liquids and gels.

medium blue, 100 Pcs.

156 102 1 Pckg. € X



SIMPLEE®
MICRO
APPLICATOR

Flocked micro applicators for precision
application of liquids and gels,
bendable. Length: 10 cm.

Colour 100 Pcs.


blue (regular) **156 105**

yellow (fine) **156 106**

white (x-fine) **156 107**

1 Pckg. € X





NiTi ring with tubes



156122



156123



156124

SIMPLEE® SECTION MATRICES

Sectional matrix system to create perfect contact points for Class II fillings, as well as to anatomically accurately form a tooth shape in the molar area.

Special silicone tubes on the ends of the ring provide a secure attachment to the matrix, while ensuring perfect matrix adaptation between the teeth.

Nickel-titanium reinforced ring with plastic grip details.

Ring and silicone tubes can be autoclaved up to 134°C.

Matrices are made of stainless steel.

Art	Ref
Kit	156 120
2 NiTi rings, 6 silicone tubes, 40 molar matrices, 40 pcs premolar och 20 pcs molar+extension	
Ring kit	156 121
2 NiTi rings, 6 silicone tubes	
Molar + extension	156 122
25 pcs matrix+extension, 0,04mm	
Molar	156 123
50 pcs molar matrices, 0,04mm	
Premolar	156 124
50 pcs premolar matrices, 0,04mm	



SIMPLEE® HELIX TEST

The Simplee Helix Test is an indicator and a batch control system for testing autoclave sterilization processes.

Intended for class B and S steam autoclaves.

PCD in accordance with NEN-EN 867-5.

The indicators meet the requirements of ISO 11140-1 class 2.

Chamber and 250 test strips.

156 715



NEW

 **SIMPLEE®**
COLDSPRAY

Cooling spray to quickly freeze foam pellets,
for vitality testing and for cooling of
impression material.

CFC-free.

Outlet temperature -50°C, mint aroma.

200 ml

156 210 1 Bottle

MINT AROMA

-50 °C



SIMPLEE®
All you need. No more. But no less.